

# 2012 Myocardial Perfusion Imaging Reimbursement

2012 CPT/HCPCS	PHYSICIAN		FACILITY	
	Reimbursement Component	MPFS Amount	APC	Hospital Outpatient Payment
CPT 78452 Myocardial Perfusion Imaging SPECT, multiple studies with wall motion and ejection fraction	Professional (26-Modifier)	\$76.45	0377	\$502.51
	Technical (Tc-Modifier)	Carrier Priced		
	Global	Carrier Priced		
CPT 78492 Myocardial Perfusion Imaging PET, multiple studies at rest & stress <small>reimbursement information taken from: <a href="http://asn.org/content_12322.cfm">http://asn.org/content_12322.cfm</a></small>	Professional (26-Modifier)	\$92.76	0308	\$1,038.36
	Technical (Tc-Modifier)	Carrier Priced		
	Global	Carrier Priced		

## Attrius®

The Attrius was developed and is optimized for myocardial perfusion imaging. The Attrius has several features making it the scanner of choice for nuclear cardiologists who value high quality PET imagery and cardiovascular specific interpretation tools, assisting in accurately assessing the patient's cardiovascular disease burden.

## Key Features

- // Low cost of ownership
- // Stress/rest protocol of 30 minutes
- // Compact footprint - fits in a 15' x 20' area
- // Proprietary cardiac software package including the ability to monitor therapy
- // Large list mode buffer allowing for concurrent acquisition of flow, perfusion and function in one scan
- // Coronary artery overlay display
- // 440 pound table limit
- // Integrated INTEL chipset with reconstruction times of 7 seconds per image set
- // Open architecture for new protocol development and customization
- // Quadrax control panels allow multi directional scanner access and control

## Disclaimer:

Positron Corporation cannot guarantee coverage and/or payment for products or procedures at any particular level. Such determinations are made based on individual patient conditions and can vary widely depending on the insurer's policies. For more specific information please contact your Medicare Contractor and/or the patient's insurer.



Positron



“Since 1995, cardiac PET has been approved as a first line test in symptomatic patients.”

## contact

530 oakmont lane  
westmont, il 60559

p 317.576.0183

f 317.576.0358

attrius@positron.com

www.positron.com

256-00030 Rev 02